

## ***S.pombe* – Simple Protein Extraction**

From: Hyun Soo Kim, Keogh lab, AECOM

**NOTE 1:** As used for confirming TAP-tag integration – this protocol is NOT recommended if you want to look at chromatin components (particularly things like modified histones).

**NOTE 2:** All buffers are pre-chilled on ice unless stated otherwise.

**NOTE 3:** Major advantage of this protocol - all manipulations can be done in a single tube, increasing throughput.

1. Grow 10 ml culture to  $OD_{600} \approx 0.5$  (or even easier, collect cells from a freshly sectoried plate). Harvest by centrifugation (pellet volume  $\approx 50\mu\text{l}$ ). Remove supernatant, resuspend cells in 1 ml ddH<sub>2</sub>O (see **NOTE 2**) and transfer to eppendorf tube. Flash spin (5 sec, 12K RPM) to collect pellet and discard supernatant.
2. Add 200 $\mu\text{l}$  Modified TEG (40mM Tris-HCl pH 7.5, 1 mM EDTA, 10% Glycerol, 0.1% NP-40, 150 mM NaCl, 1  $\mu\text{g}/\text{mL}$  leupeptin, 1 $\mu\text{g}/\text{mL}$  aprotinin, 1 $\mu\text{g}/\text{mL}$  antipain, 1 $\mu\text{g}/\text{mL}$  pepstatin A, 1 mM PMSF (see **NOTES 4 & 5**) and resuspend by pipetting. Add acid-washed glass beads to level of the meniscus (about 200-250 $\mu\text{l}$ ).

**NOTE 4:** Leupeptin, aprotinin, antipain, pepstatin A stocks all 1000 x (1mg / ml) at -20°C.

**NOTE 5:** PMSF is highly toxic (**LD50** <500mg / kg) - take care when weighing it out (stock 100mM in MeOH). Add to working solutions just before use: unstable in aqueous: half-life 110 min at pH 7, 35 min at pH 8.

3. Vortex 5 - 6 cycles (1 min vortex, 1 min on ice;  $\approx 90\%$  cell breakage).
5. Microfuge 5 sec / 1K rpm (low speed spin to collect glass beads). Add 200  $\mu\text{l}$  2X reducing loading buffer (see below) and vortex 10 sec. Microfuge 5 sec / 1K rpm. Boil 5 min 95°C. Spin 12K / 5min / RT°C. Resolve 20-30  $\mu\text{l}$  on an SDS-PAGE mini gel of the appropriate % resolving gel. **It is not necessary to transfer the S/N to a new eppendorff tube – just remove the supernatant above the glass beads / junk.**

### **2 x reducing loading buffer -**

60mM Tris pH 6.8  
2% SDS  
10% Glycerol  
0.2% Bromophenol Blue  
100mM DTT