

## One step TAP-Tag purification

**NOTE 1:** This is the one step protocol which gives great yields and has been used to purify enzyme complexes for *in vitro* functional analyses. If you want super-pure samples for MS-studies use the two step protocol (which has a much lower yield, but higher purity). I've used this protocol to purify kinases, phosphatases and HATs for *in vitro* analyses: all is the same until the dialysis buffer of choice.

**Modification of protocol used in:** Keogh *et al* (2006) *Genes Dev* 20:660-5.

### Making extracts

1. Grow 1L cells in appropriate medium to  $OD_{\lambda 600} < 1.0$  (exponential phase, etc). All steps below are with pre-chilled solutions on ice unless stated otherwise.
2. Collect cells by centrifugation. Successively wash with 100mls ddH<sub>2</sub>O and 50mls Lysis buffer (LB<sup>o</sup>, + inhibitors; LB<sup>+</sup>). Transfer cells to a 50ml Falcon and collect by centrifugation. Estimate the wet weight of the pellet ( $\approx 5\text{g/L}$  for 1L culture  $OD_{\lambda 600} \approx 1.0$ ). Store pellets @  $-80^{\circ}\text{C}$  (stable for months).
3. Thaw pellet on the bench (RT<sup>o</sup>) and then place on ice. Add 5mls LB<sup>+</sup> and resuspend completely. Add an equal volume of glass beads and place on ice for 5-10mins to chill. Lyse with 5 x vortex pulses (1 min vortex, 1 min on ice). Centrifuge 10K, 10m,  $4^{\circ}\text{C}$ . This is done in the preparative centrifuge with the Sorvall SS-34 rotor. You can spin the 50ml Falcon in the SS-34 rotor if you do not use the lid.
4. Collect the straw-colored supernatant to a Beckman untracentrifuge tube (26.3ml tube) on ice. Add another 5mls LB<sup>+</sup> to the lysate tube and resuspend completely with another 1 min vortex. Centrifuge 10K, 10m,  $4^{\circ}\text{C}$  and pool the supernatant in the ultra tube. Balance the ultra tubes and centrifuge 45K, 30m,  $4^{\circ}\text{C}$ . Collect the supernatant to a 15ml Falcon avoiding the crap at the bottom and the lipid payer at the top.
5. Using the Biorad Dc assay, estimate the concentration of the extract (usually  $\approx 10\text{mgs/ml}$ ). Given that you can collect  $\approx 8\text{mls}$  of extract this gives you  $\approx 80\text{mgs}$ (ish). Note the concentration, make 1ml aliquots and store @  $-80^{\circ}\text{C}$ .

### Purification

6. Thaw  $\approx 12\text{mgs}$  WCE on ice. Make to  $4\text{mg/ml}$  with LB<sup>+</sup> and dispense into 3 x 1.5ml ependorff tubes. Add 25 $\mu\text{l}$  IgG Agarose (Sigma, equilibrated in LB<sup>o</sup>) per tube and incubate overnight /  $4^{\circ}\text{C}$  with rotation.
7. Collect beads by centrifugation (1.5K, 2min, benchtop). Wash with 3 x 1ml LB<sup>+</sup> (NO PMSF) and pool all beads from the three ependorff tubes to one tube with 1 ml TEV buffer (TEV b<sup>o</sup>). Collect by centrifugation and remove as much supernatant as possible. Make to

final volume  $\approx 200\mu\text{l}$  with TEV b<sup>0</sup>. Remove  $\approx 5\mu\text{l}$  of the slurry to a new tube to confirm efficient digestion (before).

8. Add  $8\mu\text{l}$  Jess-made rTEV (see **NOTE 2**). **NB.** This volume of rTEV is accurate for the batch available @ 5/28/05; see pH13.56. The activity of new batches will have to be determined empirically. Digest  $4^{\circ}\text{C}$ , 4h on rotator. Remove  $\approx 5\mu\text{l}$  of the slurry to a new tube to confirm efficient digestion (after). Collect the supernatant through a Promega Wizard column. Wash beads with  $40\mu\text{l}$  TEV b<sup>0</sup> and pool with supernatant.

**NOTE 2:** TEV, Tobacco Etch Virus NIa protease. Super expensive and available from Invitrogen (Cat# 10127-017; 1000U;  $10\text{U}/\mu\text{l}$ ). Protocol for making rTEV NIa protease available at <http://mckeogh.googlepages.com/protocols>. Each batch of rTEV should be titrated for activity. Below % substrate hydrolysis data from Invitrogen.

Time	4°C	16°C	21°C	30°C
0.5h	34	58	56	77
1h	58	80	78	90
2h	71	99	99	99
3h	84	99	99	99

9. Transfer the supernatant to 2 x Slide-a-Lyzer mini dialysis units (3500MW cut-off, Pierce #69550; see **NOTE 3**). Dialyze for 90m,  $4^{\circ}\text{C}$  against  $\approx 400\text{ml}$  pre-chilled buffer of choice (sample buffers for PPases, HATs and HDACs are given below). Collect dialysate and store in aliquots at  $-80^{\circ}\text{C}$ .

**NOTE 3:** Be careful to avoid piercing the dialysis membrane with the pipette. Ensure there's nothing suspicious about the volume in the dialysis unit (I've seen a significant number leak). To check for leaks add  $100\mu\text{l}$  buffer to the unit and place it in a microfuge tube. Spin 5sec 1000rpm and discard those that allow buffer through.

## Enzymatic reactions

In each case I've tested 1, 3 and  $5\mu\text{l}$  of purified complexes and shown dose-dependence. Be sure to use an untagged strain (usually YF336 / BY4741) purified in the same manner as a control. The form the template takes is obviously up to you.

**HATs:** (example from <sup>1</sup>, see MCK lab notebook pH13.57). *In vitro* HAT reactions were performed for one hour at  $30^{\circ}\text{C}$  ( $25\mu\text{l}$  reactions containing 100ng histone substrates and  $2\mu\text{M}$  Acetyl CoA (Na salt, A2056 (mw 809.6); 1mg (\$15) to  $1125\mu\text{l} = 1\text{mM}$ ); conditions derived from <sup>2</sup>), resolved by 15% SDS-PAGE, and specific acetylation of Htz1-K14<sup>Ac</sup> determined by immunoreactivity with a specific antibody.

**PPAses:** (example from <sup>3</sup>). Phosphatase reactions with  $\gamma$ H2AX/H2B substrate were performed for one hour at 30°C (conditions derived from <sup>4</sup>), resolved by 15% SDS-PAGE, and relative phosphatase activity determined by immunoreactivity with anti- $\gamma$ H2AX.

**HDACs:** (#1, example from <sup>5</sup>; Sin3/Rpd3). *In vitro* HAT reactions were performed for one to three hours at 30°C. 25 $\mu$ l reactions contain 100ng acetylated histone substrates and HDAC under test. Samples are resolved by 15% SDS-PAGE, and specific acetylation of Htz1-K14<sup>Ac</sup> determined by immunoreactivity with a specific antibody.

**HDACs:** (#2, example from <sup>6</sup>; Sin3/Rpd3). *In vitro* HAT reactions were performed for up to one hour at 37°C (do a kinetic time course: 5, 10, 20, 40, 60m). 25 $\mu$ l reactions contain 100ng acetylated histone substrates and HDAC under test. Samples are resolved by 15% SDS-PAGE, and specific acetylation of Htz1-K14<sup>Ac</sup> determined by immunoreactivity with a specific antibody.

## Buffers

LB°     20mM Tris.Cl pH 7.6  
         10% Glycerol  
         1mM EDTA  
         200mM KoAc  
         1mM DTT  
         Protease inhibitors as STD  
         Phosphatase inhibitors if desired (1M NaF = 1000x; 0.5M Na<sub>3</sub>VO<sub>4</sub> = 1000x, 1M β-Glycerophosphate = 200X)

TEV b°     50mM Tris.HCL pH 8.0  
         1mM DTT  
         0.5mM EDTA  
         Protease inhibitors as STD (but DO NOT ADD PMSF)

HAT b°     50mM Tris.HCL pH 8.0  
         10% Glycerol  
         10mM butyric acid  
         1mM DTT  
         1mM PMSF (+ other protease inhibitors)  
         HAT b° derived from <sup>2</sup>, see MCK notebook pH13.57.

PPase b°     20mM Tris.HCL pH 7.0  
         0.1mM EDTA  
         1mM DTT  
         0.01% Brij-35  
         1mM PMSF (+ other protease inhibitors)  
         PB1 derived from <sup>4</sup>.

HDAC b° 1     20mM Tris.HCL pH 7.65  
         100mM NaCl  
         10% Glycerol  
         1mM PMSF (+ other protease inhibitors)  
         Derived from <sup>5</sup> against Sin3/Rpd3.

HDAC b° 2     10mM HEPES-KOH pH 7.8  
         75mM Kcl  
         0.1% NP40  
         5% Glycerol  
         2mM MgCl<sub>2</sub>  
         5mm DTT  
         0.25mg/ml BSA  
         1mM PMSF (+ other protease inhibitors)  
         Derived from <sup>6</sup> against Sin3/Rpd3 (I prefer the look of this one).

2xSDS-PAGE b°      120mM Tris pH 6.8  
                             2% SDS  
                             20% Glycerol  
                             0.2% bromophenol blue  
                             200mM DTT

## References

1. Keogh, M.-C., Mennella, T. A., Sawa, C., Berthelet, S. *et al.* The *Saccharomyces cerevisiae* histone H2A variant Htz1 is acetylated by NuA4. *Genes Dev* **20**:660-5 (2006).
2. Mizzen, C. A., Brownell, J. E., Cook, R. G. & Allis, C. D. in *Methods in Enzymology: Chromatin* (eds. Wasserman, P. M. & Wolffe, A. P.) (Academic Press, 1999).
3. Keogh, M.-C., Kim, J.-A., Downey, M., Fillingham, J. *et al.* DNA damage checkpoint recovery requires a phosphatase complex acting on  $\gamma$ H2AX. *Nature* **439**:497-501 (2006).
4. Zabrocki, P., Swiatek, W., Sugajska, E., Trevelein, J. M. *et al.* The *Saccharomyces cerevisiae* type 2A protein phosphatase Pph22p is biochemically different from mammalian PP2A. *Eur J Biochem* **269**, 3372-82 (2002).
5. Kadosh, D. & Struhl, K. Histone deacetylase activity of Rpd3 is important for transcriptional repression *in vivo*. *Genes Dev* **12**, 797-805 (1998).
6. Vermeulen, M. & Stunnenberg, H. G. An *in vitro* assay to study the recruitment and substrate specificity of chromatin modifying enzymes. *Biol Proced Online* **6**, 157-62 (2004).

## Modification of protocol used in:

Keogh *et al* (2006) *Genes Dev* **20**:660-5.

### *In vitro* HAT assays

Whole cell extracts (WCEs; untagged, Esa1.TAP, Gcn5.TAP or Elp3.TAP) were prepared as described (Keogh *et al.* 2003). Each HAT complex was affinity purified from 10mg WCE with IgG agarose (Sigma) (5 $\mu$ l resin / mg WCE). After overnight incubation at 4°C the bead complexes were extensively washed with Lysis Buffer (LB: 20mM Tris.Cl pH 7.6, 10% Glycerol, 200mM KoAc, 1mM EDTA, 1mM DTT + protease inhibitors), resuspended to 200 $\mu$ l final in TEV buffer (TB: 50mM Tris.HCL pH 8.0, 1mM DTT, 0.5mM EDTA) and complexes cleaved from the beads with recombinant TEV protease (4 hours, 4°C). The supernatant was collected, dialyzed for one hour against HAT buffer (HB: 50mM Tris.HCL pH 8.0, 10% Glycerol, 10mM butyric acid, 1mM DTT, 1mM PMSF) and stored in aliquots at -80°C. *In vitro* HAT reactions were performed for one hour at 30°C (25 $\mu$ l reactions containing 100ng histone substrates and 2 $\mu$ M Acetyl CoA; conditions derived from (Mizzen *et al.* 1999)), resolved by 15% SDS-PAGE, and specific acetylation of Htz1-K14 determined by immunoreactivity with anti-Htz1 K14<sup>Ac</sup>.