

Recombinant Ppr2 (TFIIS)

For *in vitro* Transcription (IVT) assays with RNAPII

NOTE 1: For general notes on recombinant protein expression see the *E.coli* lab protocols section.

NOTE 2: The eukaryotic RNAPII transcription elongation factor TFIIS is encoded by non-essential **PPR2** in *Saccharomyces cerevisiae*. *ppr2Δ* strains have normal growth on YPD but are hypersensitive to 6AU (>10μg/ml) and MPA (15μg/ml) (see: **Media Additives** protocol). In RNAPII IVTs *ppr2Δ* extracts are defective but can be easily repaired to WT levels by addition of recombinant protein: thus a great positive control for these assays. TFIIS stimulates the intrinsic 3'-5' nuclease activity of RNAPII, allowing the removal of mis-incorporated nucleotides during chain extension (Thomas *et al* (1998) *Cell* **93**:627).

See: Mkeogh notebook pH8.30 – recombinant production
Mkeogh notebook pH8.44 – Use in RNAPII IVTs

1. Transform expression plasmid pET15b-Ppr2 (BE87 / b745; **Amp^R**) into *E.coli* strain BL21/DE3.
2. Grow 100ml culture overnight at 37°C (LB^{Amp}; LBA). Sub-innoculate next morning 1/100 into 2L LBA and grow to OD₆₀₀ ≈ 0.6 at 37°C. Chill the cultures 15 minutes on ice and induce with 1mM IPTG at RT° for three hours.

NOTE 3: All steps from this point are performed at 4°C, and all buffers contain protease inhibitors as standard.

3. Harvest cells by centrifugation, wash pellet with ice-cold ddH₂O and transfer to a 50ml Falcon tubes in (can freeze pellets -80°C at this point). Resuspend in 30mls **buffer A** (4°C) and sonicate on ice (3 x 20s pulses 50% output, 40 duty cycle). Add Triton X-100 to 1% and mix by inversion.
4. Clarify by centrifugation (10K, 15m, 4°C) and transfer the supernatant to a 50ml falcon tube on ice. Add 2ml Ni-NTA resin (Qiagen) (pre-equilibrate in **buffer A**) and incubate at 4°C overnight with gentle rolling.
5. Collect beads by centrifugation (1.5K, 5m, 4°C) and wash in batch with 20mls **buffer A** and 20mls **buffer E / 20mM Imidazole**. Collect after each resuspension by centrifugation and discard the supernatants (keep some to test by Western analysis). Transfer the beads to a disposable column and elute with **buffer E / 200mM Imidazole**. Collect 500μl fractions.

6. Analyze fractions by Biorad Dc and SDS-PAGE (**expected size $\approx 35\text{kD}$**). Estimate the cleanliness of the prep (how many bands visible on commassie staining SDS-PAGE) and yield per L (**NB.** Original prep (mkeogh notebook pH8.30) $>90\%$ pure).

7.1 (Option 1) Pool positive fractions and dialyze against 2 x 1L *Transcription buffer B* (**Tx-b^oB**; 2 hrs per, 4°C). Compare concentration (by Biorad Dc) before and after dialysis (to ensure no holes in the membrane) and store in small aliquots at -80°C (recommend $500\mu\text{g/ml}$).

7.2 (Option 2) Positive fractions were pooled and the buffer changed to *Transcription buffer B* (**Tx-b^oB**) by three successive 10ml concentrations in a Centriprep YM-10. The final concentrate was made to $500\mu\text{g/ml}$ in **Tx b^oB** and stored in aliquots at -80°C .

Buffers:

Buffer A -

10mM Tris.HCl (pH 8.0)
10% Glycerol
500mM NaCl
0.1% Tween 20
10mM Imidazole
10mM β -ME
Protease inhibitors as standard (incl. 1mM PMSF)

Buffer E -

20mM HEPES.KOH (pH 7.6)
20% Glycerol
100mM NaCl
1mM EDTA
1mM DTT
Protease inhibitors as standard (incl. 1mM PMSF)
Imidazole as indicated

Transcription Buffer B -

20mM HEPES.KOH (pH 7.5)
20% Glycerol
150mM KAc
10mM MgAc
10mM EGTA
5mM DTT
Protease inhibitors as standard (incl. 1mM PMSF)

For a manuscript methods section:

Recombinant Ppr2 (TFIIS): E. coli strain BL21/DE3 was transformed with pET15b-Ppr2. A 1L culture was grown at 37°C to $\text{OD}_{600} \approx 0.5$ and induced by the addition of 1mM IPTG. Cells were further cultured for three hours at room temperature. All steps from this point

on were performed at 4°C. Lysates were prepared and Ni-NTA purified as above. Positive fractions were pooled and the buffer changed to Transcription buffer B (Tx b°B) (see *In vitro transcription* protocol) by three successive 10ml concentrations in a Centriprep YM-10. The final concentrate was made to 500µg/ml in Tx b°B and stored in aliquots at -80°C.