

Lys counter-selection with α AAA

(α -AminoAdipic Acid (α AAA): US Biological, A1374-09: 5g \$104)

NOTE 1: Many yeast labs use 5-FOA as a counter-selection agent to URA⁺, but few use α AAA against LYS⁺, usually because they've heard it doesn't work. Not so: it does, just not as well as FOA (*Anindita Basak, Personal Communication*). However when URA is unavailable, LYS is used in our lab.

NOTE 2: *S.cerevisiae* LYS2 encodes α -Aminoadipate reductase, LYS5 a phospho-pantetheinyl transferase that activates Lys2p (Chattoo *et al*, 1979a). Individual *lys2* Δ or *lys5* Δ strains are viable, Lys⁻ and resistant to α -AAA. Contrast with Lys⁺ strains which convert α AAA to a toxic intermediate.

NOTE 3: When Lys⁻ strains are available they're usually *lys2* Δ (ORF = 4179bp). The pRS317 yeast shuttle vector (MKF118) complements LYS2 but the empty vector is already 8.3kb, complicating downstream cloning. The LYS5 ORF is much smaller (819bp) which facilitates cassette construction for PCR (eg: *geneX* Δ ::Lys5.MX) (Ito-Harashima & McCusker, 2004)

1. α -AminoAdipic acid (α AAA) stock (4% in ddH₂O: will be used as 0.2% final)
Stock = 2g / 50ml (4%) in ddH₂O. Adjust pH to 6.0 with 5M KOH (see **NOTE 4**)
Filter sterilize (0.2 μ) and store at 4°C. Use @ 0.2% final (50 ml stock / L)

NOTE 4: The pH-ing to 6.0 step is absolutely required to get the α AAA into solution. It may help to initially overshoot (to \approx 8.0) and as more α AAA dissolves come closer to the desired 6.0. This will take hours: plan accordingly.

2. α AA media (per L: see **NOTE 5** – optionals indicated by a *)

YNB (No NH ₄ SO ₄)	3.0 gm
Inositol (200 mM stock)	2.0 ml
Lys (100mM stock)	10.0 ml
* Ura (20mM stock)	10.0 ml
* His (100mM stock)	3.0 ml
* Trp (40mM stock)	10.0 ml
* Leu (100mM stock)	16.7 ml
* Met (100 mM stock)	14.76 ml
Agar	20.0 gm
Glucose	20.0 gm
ddH ₂ O	to 950 ml

Autoclave w/o α AAA. Cool to \approx 50°C and add 50ml 0.2 μ filter sterilized 4% α AAA / L

NOTE 5: On this media the α AAA will be used as a principle N-source by the yeast (hence the reason to use NH₄SO₄-free YNB). In LYS⁺ strains this leads to the generation of a toxic intermediate of the lysine biosynthetic pathway while *lys2* or *lys5* mutants grow unhindered. Be careful NOT to include any additional N- sources to the media: this includes the lab -6 or -4 amino acid powders (**which contain**

glutamate) or YNB with NH_4SO_4 . As above use only those specific amino acids that supplement the auxotrophies of the strains used in the experiment – although not tested it is likely you can also supplement with Adenine.

3. α AAA is used as a counter selection agent (eg. to select for loss of a plasmid encoding a Lys2 or Lys5 locus in an otherwise *lys2* Δ or *lys5* Δ strain). This is preferably done by **Spotting**: follow the standard lab protocol (<http://mckeogh.googlepages.com/protocols>).

NOTE 6: It is ESSENTIAL to replica plate: LYS^+ strains will grow through on the first round of selection on an α AAA plate and can in fact be rescued if you replica plate back onto YPD. Two rounds of selection on α AAA media eliminates the majority of background – we have tried additional rounds but it doesn't help: those colonies that grew through on the second plate from an initial LYS^+ parent were now stable on α AAA, presumably because they had acquired a mutation at *LYS2* or *LYS5*.

NOTE 7: Those familiar with the FOA / URA system will not be overly impressed with α AAA / LYS counter-selection: the requirement for replica plating and relatively high backgrounds are significant issues. However it is a useful backup option if URA is otherwise unavailable.

REFERENCES

Chattoo BB *et al* (1979a) Selection of *lys2* mutants of the yeast *Saccharomyces cerevisiae* by the utilization of α -Aminoadipate. *Genetics* **93**: 51-65

Chattoo BB *et al* (1979b) Patterns of genetic and phenotypic suppression of *lys2* mutations in the yeast *Saccharomyces cerevisiae*. *Genetics* **93**: 67-79

Ito-Harashima S & McCusker, JH (2004) Positive and negative selection *LYS5MX* gene replacement cassettes for use in *Saccharomyces cerevisiae*. *Yeast* **21**: 53-61

Sikorski RS & Boeke JD (1991) *In vitro* mutagenesis and plasmid shuffling: from cloned gene to mutant yeast. *Methods in Enzymology* **194**: 302-318

Zaret KS & Sherman F (1985) α -Aminoadipate as a primary nitrogen Source for *Saccharomyces cerevisiae* mutants. *J Bacteriology* **162**: 579-83