

## **Purification of His-TEV Protease**

**NOTE 1:** TEV from Invitrogen (or whatever they're called these days) is super expensive (Cat# 10127-017; 1000U; 10U/ $\mu$ l). Following is protocol from Jessica Dermody. Plasmid is pTPSN from Jennifer Doudna (Lucast *et al* 2001 *BioTechniques* **30**:544) ([Buratowski lab](#) F587/b2664. [Keogh Lab](#) MKF308/mb323). Expresses HIS-tagged TEV (Tobacco Etch Virus) NIa protease with S219N mutation that removes a self cleavage site. Most expressed protein is insoluble. See reference for construction and expression steps.

**Induction:** rTEV-HIS6 under regulation of a T7 Promoter (use: *E.coli* Strain: BL21/DE3 *pLysS* or suchlike). Standard growth / induction (@ OD<sub>600</sub> 0.4, add 0.1mM IPTG / 4h, RT°C). Harvest and store at -80°C until ready to proceed.

1. Resuspend induced *E.coli* pellets in 5x volumes Ni-NTA lysis buffer  
50mM Tris-HCl, pH 8.0  
150mM NaCl  
10mM Imidazole, pH 8.0  
10mM  $\beta$ ME  
150 $\mu$ M PMSF
2. Sonicate to lyse cells: 3 x 30 sec on, 30 sec off, power 4  
(don't try extra sonication pulses to increase yield – it gives higher backgrounds)
3. Spin down lysate for 30min at 15000rpm at 4°C in SS34 rotor.
4. Mix clear lysate with prewashed Ni-NTA resin (1-2mL of slurry per liter of culture).
5. Incubate lysate plus resin on rocker at 4°C for 30-60min.
6. Spin down lysate+resin for 5min at 1000rpm and remove supernatant and save as flowthrough. Pour resin into small Biorad column. **DO ALL THIS IN THE COLD ROOM!!!**
7. Wash resin 5x column volume (cv) with wash buffer, and collect in 1 cv fractions.  
50mM Tris-HCl, pH 8.0  
150mM NaCl  
20mM Imidazole, pH 8.0  
10mM  $\beta$ ME  
150 $\mu$ M PMSF
8. Elute 5x column volume with elution buffer, and collect in 1 cv fractions.  
50mM Tris-HCl, pH 8.0  
150mM NaCl  
250mM Imidazole, pH 8.0  
10mM  $\beta$ ME  
150 $\mu$ M PMSF

9. Run 15% SDS-PAGE to determine which elution fractions contain the protease
10. Dialyze protease into:  
20% glycerol  
50mM Tris-HCl, pH 8.0  
1mM EDTA, pH 8.0  
5mM DTT
11. Store in 20µl aliquots at –80°C

**NOTE 2:** You will have to determine the home-made TEV amount empirically (see MCK lab notebook pH13.56). It is handy to have some of the invitrogen stuff to hand when making the first batch to see how well you're doing. The enzyme is usually so active that four hours at 4°C should go to completion. It is possible to determine this by  $\alpha$ PAP (Sigma) reactivity in westerns – you're looking for the disappearance of the TAP-band. You can also use the  $\alpha$ CBP antibody (although I'm not a huge fan of the latter). Below % substrate hydrolysis data from Invitrogen.

Time	4°C	16°C	21°C	30°C
0.5h	34	58	56	77
1h	58	80	78	90
2h	71	99	99	99
3h	84	99	99	99

**rTEV purified with this protocol was used in:**

Keogh *et al* (2006) *Nature* **439**:497

Keogh *et al* (2006) *Genes Dev* **20**:660-5