

## Media Additives

### *Antibiotics, counterselection agents & inhibitors*

**NOTE 1:** Use in conjunction with the **Media Recipes** protocol.

**NOTE 2:** **NEVER** autoclave media with an additive. Generally an additive should not be added to agar-based media until cooled to  $\approx 50^{\circ}\text{C}$  (just about cool enough to hold in your palm). Mix well (solubility is an issue with some of these things), 70% EtOH spritz and flame to remove bubbles (**careful**) and pour plates (generally a little thinner than normal to keep costs down).

### Antibiotics (*E.coli*)

Ampicillin: (Amp) 1000 x stock = 50mg/ml in ddH<sub>2</sub>O (>99% of the time – some older plasmids have special requirements: **stringency issue**). 0.2 $\mu$  filter sterilize and store 1ml aliquots at  $-20^{\circ}\text{C}$ . Plates stable for >6 months at  $4^{\circ}\text{C}$ . Mixed media can be used after a week or two at RT $^{\circ}\text{C}$  (but storage at  $4^{\circ}\text{C}$  is recommended).

**NOTE:** Amp is bacteriostatic rather than bactericidal.  $\beta$ -lactamase (the Amp<sup>R</sup> product) is secreted into the media around Amp<sup>R</sup> colonies, so satellites can develop if the plates are incubated for too long.

Kanamycin: (Kan) 1000 x stock = 50mg/ml in ddH<sub>2</sub>O (>99% of the time – some older plasmids have special requirements: **stringency issue**). 0.2 $\mu$  filter sterilize and store 1ml aliquots at  $-20^{\circ}\text{C}$ . Plates stable for >6 months at  $4^{\circ}\text{C}$ . Mixed media can be used after a week or two at RT $^{\circ}\text{C}$  (but storage at  $4^{\circ}\text{C}$  is recommended).

**NOTE:** Kan is bactericidal – satellites not a problem

Chloramphenicol: (Chl) 1000 x stock = 34mg/ml in EtOH. Store at  $-20^{\circ}\text{C}$ . Usually used to drive-up the yield of some low-copy plasmids: rarely done nowadays.

Tetracycline: (Tet) 1000 x stock = 5mg/ml in EtOH. Store at  $-20^{\circ}\text{C}$ . Usually used to maintain the F' plasmid of some color-selectable strains (check the genotype maps) prior to making competent cells.

### Antibiotics (*S.cerevisiae* / *S.pombe*)

G418: (G418, geneticin sulfate, neomycin sulfate; eg. American Bioanalytical AB05057) 250 x stock = 50mg/ml in ddH<sub>2</sub>O. 0.2 $\mu$  filter sterilize and store 1ml aliquots at  $-20^{\circ}\text{C}$  for lab stocks. Plates (usually @ 200 $\mu\text{g/ml}$  in YPD) stable for >6 months at  $4^{\circ}\text{C}$ .

**NOTE:** When transforming the KanMX cassette into yeast it is strongly advised to replica plate 24hrs after transformation. If this is omitted there is usually a high background and it can be difficult to identify true G418<sup>R</sup> colonies.

Noursethricin: (clonNAT, Nat) 1000 x stock = 100mg/ml in ddH<sub>2</sub>O. 0.2μ filter sterilize and store 1ml aliquots at -20°C. Plates (usually made in YPD) stable for >6 months at 4°C.

**NOTE:** (from Werner Bioagents, <http://www.webioage.com/seite5.html>) ClonNAT, trade name for nourseothricin (a complex of the streptothricins F and D), produced by *Streptomyces noursei*, as a dihydrogen sulphate. Inhibits bacteria, mycobacteria, mycoplasmas, yeasts, viruses and plants by inhibiting ribosomal protein synthesis and inducing miscoding.

## Other additives (*S.cerevisiae*)

α-aminoadipic acid: (DL-αAAA; eg. US Biological A1374-09 5g \$104) Use at 0.2% final.

Used to counterselect *LYS*2<sup>+</sup>.

Stock = 2g / 50ml (4%) in ddH<sub>2</sub>O. Adjust pH to 6.0 with 5M KOH, 0.2μ filter sterilize and store at 4°C. Use 50ml of this per L media (0.2% final): **see specific protocol on LYS Counterselection (<http://mckeogh.googlepages.com>)**

**NB. Must add *LYS*<sup>+</sup> to media (3.3ml 100mM stock / L)**

**NOTE:** Solubility is an issue and the pH-ing to 6.0 step is absolutely required to get the αAAA into solution. It may help to initially overshoot (to ≈ 8.0) and as more αAAA dissolves come closer to the desired 6.0. This will take hours: plan accordingly. 5-FOA is more commonly used as a counterselection agent

**CLASS:** Counterselection agent

5-FOA: (5-FluorO-Orotic acid; eg. US Biological F5050 10g \$199) Use at 1mg/ml. Used to counterselect *URA*3<sup>+</sup>. **Use *URA*<sup>+</sup> synthetic complete plates** (the conversion of 5-FOA to the toxic metabolite requires a functional *URA*3, and you're selecting against this) to shuffle out a *URA* plasmid.

**NOTE:** Solubility is an issue. Add appropriate amount for final volume to 2 x synthetic complete medium (+ URA). Place on shaker for an hour or two at 37°C. Add an equal volume of >70°C 2 x Agar (40g/L in ddH<sub>2</sub>O), Mix well, 70% EtOH spritz and flame to remove bubbles (**careful**), and pour plates (generally a little thinner than normal to keep costs down). Plates stable for >6 months at 4°C.

**CLASS:** Counterselection agent

6AU: (6-AzaUracil) Stock = 100mg/ml in DMSO. Store at -20°C. Commonly used at 75μg/ml final in *URA*<sup>-</sup> synthetic complete plates. Plates stable for >6 months at 4°C.

Transcription elongation inhibitor; works by depleting nucleotide pools (see also MPA). The massive amount of Uracil in *URA*<sup>+</sup> plates kills this. Strain should be *URA*<sup>+</sup> (ie. *URA*3 in the genomic locus, or transformed with a cen/ars *URA*3 plasmid, eg. pRS316).

**NOTE:** While commonly used at 75μg/ml, it is recommended to test a variety of concentrations: 0, 10, 25, 50, 75 and 200 (all μg/ml final). *ppr2Δ* (TFIIS) shows a

sensitivity at 10, although most don't show problems until > 25 (Keogh *et al* (2004) *Mol Cell Biol* **23**:7005).

**CLASS:** Transcription elongation inhibitor

Benomyl: (Ben) Commonly used 15µg/ml final in YPD plates. **Horribly insoluble.**

**Used in:** Krogan *et al* (2004) *PNAS* **101**:15313; Keogh *et al* (2006) *Genes Dev* **20**:660.

**CLASS:** Microtubule destabilizer, chromosome stability analyses.

Bleomycin: (bleo) Commonly used at 5µg/ml final in YPD plates. **Expensive.** Generates free radicals and induces DNA lesions similar to those caused by ionizing radiation.

**Used in:** Keogh *et al* (2006) *Nature* **439**:497.

**CLASS:** DNA damaging agent / genotoxin

Camptothecin: (CPT) Commonly used at 20µM final (7mg/L, 100mg bottle) in YPD plates. Topoisomerase inhibitor that induces DSBs by causing replication forks to stall.

**CLASS:** DNA damaging agent / genotoxin

HydroxyUrea: (HU) Use at 100mM final (Yes, that is milliMolar), commonly in YPD medium or plates. Inhibitor of dNTP synthesis, leads to DNA replication fork collapse. Also used to induce a cell-cycle arrest (early S-phase).

**Used in:** Keogh *et al* (2006) *Nature* **439**:497; Keogh *et al* (2006) *Genes Dev* **20**:660.

**CLASS:** DNA damaging agent / genotoxin, cell-cycle arrest

MethaneMethyl Sulfonate: (MMS) 100% stock in chemical cabinet. Use at 0.05 - 0.1% final (commonly in YPD medium or plates). Alkylating agent.

**Used in:** Keogh *et al* (2006) *Nature* **439**:497; Keogh *et al* (2006) *Genes Dev* **20**:660.

**CLASS:** DNA damaging agent / genotoxin

MycoPhenolic acid: (MPA) Commonly used at 15µg/ml final in URA<sup>-</sup> synthetic complete plates. Transcription elongation inhibitor; works by depleting nucleotide pools (see also 6AU). The massive amount of Uracil in URA<sup>+</sup> plates kills this. Strain should be URA<sup>+</sup> (ie. *URA3* in the genomic locus, or transformed with a cen/ars *URA3* plasmid, eg. pRS316).

**NOTE:** When working with a possible transcription elongation factor, it is recommended to test for MPA sensitivity (it usually mirrors 6AU), but the latter is more commonly used for in-depth experiments (Keogh *et al* (2004) *Mol Cell Biol* **23**:7005).

**CLASS:** Transcription elongation inhibitor

Nocodazole: (noc) Make a 1.5mg/ml stock in DMSO; this is 100x. (Sigma M1404 >99%TLC; 10mg \$54, 50mg \$215). Destabilizes microtubules preventing sister chromatid separation: G<sub>2</sub>/M phase.

**NOTE:** Nocodazole escape can occur (particularly with checkpoint mutants). Arresting cells for too long is eventually lethal. Nocodazole arrest is strain dependent (too much can induce break-through, results are also problematic at 37°C).

**Used in:** Keogh *et al* (2006) *Genes Dev* **20**:660.

**CLASS:** Microtubule destabilizer, cell-cycle arrest

Sodium Butyrate: Histone DeACetylase (HDAC) inhibitor. Use at 5 – 50mM final (test at 0.5mM). Relatively cheap, but nowhere near the potency of Trichostain A (TsA).

**CLASS:** HDAC inhibitor

Trichostain A: (TsA) HDAC inhibitor (see also Sodium Butyrate). Use at 30nM – 3μM final. **Very Expensive.** Used in (Keogh *et al* (2005) *Cell* **123**:593).

**CLASS:** HDAC inhibitor