

***S. cerevisiae* acetone powder**

NOTE: Used for the pre-adsorption of polyclonal antisera, this can reduce backgrounds significantly. If you raise an Ab to protein Yfg, a powder will be made to the *YFG* deletion strain and used to adsorb a batch of the appropriate Yfg anti-sera: if you were to use a WT strain, you'd obviously just adsorb the specific signal.

1. Grow 50ml culture of the appropriate Δ strain to confluence (late log, $OD_{600} \approx 2$). Collect cells by centrifugation.
2. Resuspend cells in 1x PBS at a density of 10 – 30 OD_{600} / ml. Add zymolase ($\approx 1\text{mg}/1000$ ODs) and incubate with gentle shaking 10 – 20 min 30°C
3. Place on ice 10m. Don't worry if spheroplasting is incomplete: it's really only necessary to weaken the cell wall.
4. Flick pellet to resuspend and add 4 vol acetone (-20°C). Vortex to a single cell suspension and put on ice 30 min.
5. Collect by centrifugation and resuspend pellet in acetone (-20°C). This may take a bit of work and involve some crushing with a large wooden applicator and extreme vortexing. Put on ice 30 min.
6. Collect by centrifugation and transfer pellet to mortar. Leave until dry (should be pretty quick – acetone is very volatile but can speed up using fume hood).
7. Grind to a fine powder using pestle. Store powder at RT° in a foil-wrapped falcon tube.
8. Use powder in undiluted sera at 10mg/ml and adsorb overnight at 4°C with rotation
9. Centrifuge 12000rpm, 10min, 4°C . Remove adsorbed sera to a new ependorff and compare with pre-adsorbed sample.