

## **Salmon sperm DNA**

For use in yeast transformations and blocking agent in blotting

*Use a better grade as a carrier in nucleic acid precipitations*

1. Weigh out 400mg dried salmon sperm DNA (**eg. Sigma Sperm nuclei type II-S; S3126**)  
Cut the strands with scissors to small fibers.
2. Place DNA in a 50ml Falcon tube. Add ddH<sub>2</sub>O to 40mls (10mg/ml final). Swirl to release trapped bubbles.
3. Heat tube in 65°C bath overnight to dissolve DNA. Invert tube occasionally to aid resuspension.
4. When in solution, the DNA will be a gooey mess. Sonicate at about 50% power for a minute at a time for several cycles until the solution becomes less viscous. A good rule of thumb is when bubbles can rise easily to the top (**this is real rocket science**).
5. Pierce the lid of the Falcon with a needle and place in a beaker of boiling H<sub>2</sub>O for 15 minutes. Repeat sonication.
6. Cool and dispense 1ml aliquots. Store at -20°C.
7. To use as a carrier in yeast transformations pierce cap and boil for 5 mins. Store at 4°C for repeated usage.