

## Making electrocompetent *E.coli*

*See allied protocol for electrotransformation*

**NOTE 1:** Electrocompetent *E.coli* generally have a higher transformation efficiency (more colonies per  $\mu\text{g}$  / DNA) than chemically competent cells (**Expected  $\geq 5 \times 10^9$  cells /  $\mu\text{g}$  DNA**). On the minus side the protocol is more fiddly than heat-shocking.

**NOTE 2:** *E.coli* strains used (detailed genotypes online)

(i)	XL1-blue	Tet <sup>R</sup>	
(ii)	S28	Tet <sup>R</sup>	vv. weak pellet
(iii)	DH5 $\alpha$	Tet <sup>S</sup>	stronger pellet
(iv)	DH10b	Tet <sup>S</sup>	

1. Confirm the genotype of the strain to be made competent. Determine if it has an F' (usually required for color-selection), and thus whether the seed should be grown with Tetracycline.
2. Inoculate colony into 50ml LB (+ Tet if appropriate). Grow overnight: 37°C, 300rpm.
3. Sub 1/100 into 1L LB. Note OD<sub>600</sub>. Incubate 37°C, 300rpm. Measure OD<sub>600</sub> every 30 – 60 mins. Harvest at OD<sub>600</sub>  $\approx$  0.5 (range 0.4  $\rightarrow$  0.6). Chill culture 15 min in an ice-bath.

**NOTE 3:** ALL buffers from now on are sterile and pre-chilled on ice. Keep everything  $\leq 2^\circ\text{C}$ .

4. Transfer culture to centrifuge flasks and collect 4000rpm, 2°C, 10m.

**NOTE 4:** The quality of the *E.coli* pellet will get worse after each spin (presumably the *E.coli* contacts are ion-dependent).

5. Discard supernatant (**careful**) and resuspend pellet in 250ml ddH<sub>2</sub>O. Centrifuge 4000rpm, 2°C, 10m. Repeat this step twice.

**NOTE 5:** Pre-label ependorffs and store rack at  $-20^\circ\text{C}$ . Also get some dry-ice.

6. Resuspend pellet in 40ml 10% glycerol. Transfer cells to a 50ml Falcon tube. Centrifuge 6000rpm, 2°C, 10m. Estimate packed cell volume and add an equal volume of 10% glycerol. Dispense 200 $\mu\text{l}$  aliquots into tubes on dry-ice and store at  $-80^\circ\text{C}$ .
7. Estimate the competence (expressed as colonies per  $\mu\text{g}$  plasmid) by transforming with a known amount of plasmid. Plate numerous dilutions and count cells. Note the competence and date on the box and consider that the efficiency will drop with increasing storage, even at  $-80^\circ\text{C}$  (best to use cells within  $\leq 12$  months).

**NOTE 6:** The number of transformants increases linearly with input DNA over a wide range (5pg  $\rightarrow$  500ng). However it is important to keep the added volume of DNA as small as possible (particularly important for ligations) since arcing can occur at salt concentrations  $>1\text{mM}$ .