



Construction of **deletion** cassette with KanMX6 (Nat, Hph, Pat)



Construction of **C-terminus Tagging** cassette with KanMX6



PCR conditions for 2-step Deletion/TAP method in *S. pombe*

1. Not1 Digestion of template cassettes

pFA6a series vectors	(2~4ug)
NEBuffer 3	3ul
10X BSA	3ul
Not1	3ul
DW	to 30ul

→37°C for O/N → purification with PCR purification kit → 0.7% Gel running & EtBr Staining

2. 1st PCR (Deletion & C-terminal Tagging Cassette)

Genomic DNA (WT)	3ul	<Pombe 1 / HSK> 94/4m 94/30s 50 ^{***} /30s 72/1m (X30) 72/5m
10X buffer	5ul	
dNTP mix (2mM each)	5ul (final 0.2mM each)	
L1* (or L3, L5) (10uM)	1.5ul	
L2* (or L4, L6) (10uM)	1.5ul	
Taq (Sure-Pol ^{**})	2ul	
DW	32ul (to 50ul)	

* L1/L2, L3/L4, L5/L6 pair

** Pfu version taq polymerase

*** Try 40°C or 45°C if there is no product (usually L5/L6 product has problem because of high AT-content of L5)

→ purification with PCR purification kit → 0.7% Gel running & EtBr Staining → measure DNA concentration

3. 2nd PCR (Deletion & C-terminal Tagging Cassette)

Template cassette	50ng	<Pombe 2 / HSK> 94/4m 94/15s 40/30s 72/3m (X5) 94/15s 55/30s 72/3m (X25) 72/5m
10X buffer	5ul	
dNTP mix (2mM each)	5ul (final 0.2mM each)	
L1/L2 product (or L3/L4)	50ng	
L5/L6 product	50ng	
L1 (or L3) (10uM)	1ul	
L6 (10uM)	1ul	
Taq (Sure-Pol)	2ul	
(DMSO* for NAT cassette)	2.5ul (5%)	
DW	to 50ul	

* because of high GC-content in NAT cassette (other cassettes are fine)

→purification with PCR purification kit → 0.7% Gel running & EtBr Staining → measure DNA concentration

→Transformation to *S. pombe* cell (500ng-1ug, up to 10ul per rxn)

4. PCR screening

Genomic DNA	3ul	<Pombe 3 / HSK> 94/4m 94/30s 50/30s 72/3m (X35) 72/5m (* you can also use A1/KanMX6-5'- Rev (#505) or KanMX6-3'- For(#1134)/A3)
10X buffer	3ul	
MgCl ₂ (50mM)	1.5ul (2.5mM)	
dNTP mix (2mM each)	3ul (0.2mM each)	
A1* (or A2) (10uM)	0.6ul	
A3 (10uM)	0.6ul	
Taq (Taq-pro)	0.3ul	
DW	18ul (To 30 ul)	

→ Control sample must be included !!

Oligomers for Deletion & C-terminal tagging (forward)

spSwc6-L1 (20mer) :

spSwc6-L2 (45mer) : taattaacccggggatccgtcgacctatgagctactcaaaaattg

spSwc6-L3 (20mer) :

spSwc6-L4 (45mer) : taattaacccggggatccgtcgaccggcatatacttttcatacacc

spSwc6-L5 (46mer) : aaacgagctcgaattcatcgatgataaagttgaccttttgcttact


spSwc6-L6 (20mer) :

spSwc6-A1 (20mer) :

spSwc6-A2 (20mer) :

spSwc6-A3 (20mer) :

"add gene-specific sequences to colored bold sequences"



(L1-L2, L3-L4, L5-L6: 200~300bp)

caatTTTTgagtagctcata**atg**tttgtcaccgccattgaaca
 gttaaaaactcatcgagtat**tac**aaacagtgggggtaacttgt
 ← L2 **ccagctgcctaggggccaattaat**

aaacgagctcgaattcatcgatgata
 ← L5
 ggtgtatgaaagtatatgcc**taa**aagttgaccttttgcttact
 ccacatactttcatatacgg**att**ttcaactggaaaacgaatga
 ← L4 **ccagctgcctaggggccaattaat**

Oligomer for **Deletion** & **C-termial** tagging (**Reverse**)

spAsh2-L1 (20mer) :

spAsh2-L2 (45mer): taattaacccgggatccgtcgaccggctgattaatattgaatcc

spAsh2-L3 (20mer) :

spAsh2-L4 (45mer): taattaacccgggatccgtcgaccgtttaaatagccacgacatg

spAsh2-L5 (46mer): aaacgagctcgaattcatcgatgataacaatgatgctcttgagaa



spAsh2-L6 (20mer) :

spAsh2-A1 (20mer) :

spAsh2-A2 (20mer) :

spAsh2-A3 (20mer) :

```
"add gene-specific sequences to colored bold sequences"
```



pFA6a-KanMX6

(all pFA6a cassettes have same DAN sequences except resistant genes (Kan^R, Nat^R..)

```
1 gaacgcgggcc gccagctgaa gcttcgtacg ctgcaggctcg acggatcccc gggttaatta multi cloning site
61 aggcgcgcca gatctgttta gcttgccctcg tccccgcgcg gtcacccggc cagcgacatg
121 gagggccaga ataccctcct tgacagtctt gacgtgcgca gctcaggggc atgatgtgac
181 tgtcgcccggt acatttagcc catacatccc catgtataat catttgcatc catacatttt KanMX6-5' Rev(#505)
241 gatggccgca cggcgcggaag caaaaattac ggctcctcgc tgcagacctg cgagcagggg TEF promoter
301 aacgctcccc tcacagacgc gttgaattgt cccacgcgcg cgccctgtga gagaaatata
361 aaagggttagg atttgccact gaggttcttc ttcatatac ttcttttaa aatcttgcta
421 ggatacagtt ctacatcac atccgaacat aaacaaccat gggtaaggaa aagactcacg
481 tttcgaggcc gcgattaaat tccaacatgg atgctgattt atatgggtat aaatgggctc
541 gcgataatgt cgggcaatca ggtgcgacaa tctatcgatt gtatgggaag cccgatgcgc
601 cagagttggt tctgaaacat ggcaaaggta gcgttgccaa tgatgttaca gatgagatgg
661 tcagactaaa ctggctgacg gaatttatgc ctcttcgcac catcaagcat tttatccgta
721 ctctgatga tgcattggtta ctccactg cgatcccccg caaacagca ttccaggtat
781 tagaagaata tctgattca ggtgaaaata ttgttgatgc gctggcagtg ttctgcgcc
841 ggttgcatte gattcctggt tgtaattgtc cttttaacag cgatcgcgta tttcgtctcg KanR
901 ctcaggcgca atcacgaatg aataacggtt tgggtgatgc gagtgatttt gatgacgagc
961 gtaatggctg gcctgttgaa caagtctgga aagaaatgca taagcttttg ccattctcac
1021 cggattcagt cgtcactcat ggtgatttct cacttgataa ccttattttt gacgagggga
1081 aattaatagg ttgtattgat gttggacgag tcggaatcgc agaccgatac caggatcttg
1141 ccacccatg gaactgcctc ggtgagtttt ctcttcatt acagaaacgg ctttttcaa
1201 aatatgggtat tgataatcct gatatgaata aattgcagtt tcatttgatg ctcgatgagt
1261 tttcttaa agtactgaca ataaaaagat tcttgttttc aagaacttgt catttgata
1321 gtttttttat attgtagttg ttctatttta atcaaagtgt agcgtgattt atattttt ttt KanMX6-3' For (#1134)
1381 tcgcctcgac atcatctgccc cagatgcgaa gttaagtgcg cagaaagtaa tatcatgcgt TEF terminator
1441 caatcgtagt tgaatgctgg tcgtatact gctgtcgatt cgatactaac gccgcatcc
1501 agtttaaagc agctcgaatt catcgatgat atcagatcca ctagtggcct atgcggccgc multi cloning site
1561 ggatctgccc gtctccctat agtgagtcgt attaatttcg ataagccagg ttaacctgca
1621 ttaatgaatc ggccaacgcg cggggagagg cggtttgcgt attgggcgct cttccgcttc
1681 ctgctcact gactcgctgc gctcggtcgt tcggctgcgg cgagcgggtat cagctcactc
1741 aaaggcggtg atacggttat ccacagaatc aggggataac gcaggaaaga acatgtgagc
1801 aaaaggccag caaaaggcca ggaaccgtaa aaaggccgcg ttgctggcgt tttccatag
1861 gctccgcccc cctgacgagc atcacaaaaa tcgacgtca agtcagaggt ggcgaaaccc
```

KanMX6-5' Rev (#505): acatggggatgtatgggcta

KanMX6-3' For (#1134): tttcgcctcgacatcatct

1921 gacaggacta taaagatacc aggcgtttcc ccttggaagc tccctcgtgc gctctcctgt
 1981 tccgaccctg ccgcttaccg gatacctgtc cgcctttctc ccttcgggaa gcgtggcgct
 2041 ttctcaatgc tcacgctgta ggtatctcag ttcggtgtag gtcgttcgct ccaagctggg
 2101 ctgtgtgcac gaaccccccg ttcagcccga ccgctgcgcc ttatccggtg actatcgtct
 2161 tgagtccaac ccggttaagac acgacttata gccactggca gcagccactg gtaacaggat
 2221 tagcagagcg aggtatgtag gcggtgctac agagtctctg aagtgggtggc ctaactacgg
 2281 ctacactaga aggacagtat ttggtatctg cgctctgctg aagccagtta ccttcggaaa
 2341 aagagttggg agctcttgat ccggcaaaaa aaccaccgct ggtagcgggtg gtttttttgt
 2401 ttgcaagcag cagattacgc gcagaaaaaa aggatctcaa gaagatcctt tgatcttttc
 2461 tacgggggtct gacgctcagt ggaacgaaaa ctacagttaa gggattttgg tcatgagatt
 2521 atcaaaaagg atcttcacct agatcctttt aaattaaaaa tgaagtttta aatcaatcta
 2581 aagtatatat gagtaaactt ggtctgacag **tta**ccaatgc ttaatcagtg aggcacctat
 2641 ctacgcgac tgtctatttc gttcatccat agttgcctga ctcccgcctg ttagataaac
 2701 tacgatacgg gagggcttac catctggccc cagtgtctga atgataccgc gagaccacg
 2761 ctaccgggt ccagatttat cagcaataaa ccagccagcc ggaagggccg agcgcagaag
 2821 tggctctgca actttatccg cctccatcca gtctattaat tgttgccggg aagctagagt
 2881 aagtagttcg ccagttaata gtttgcgcaa cgttgttgcc attgctacag gcacgtgggt
 2941 gtcacgctcg tcgtttggta tggtctcatt cagctccggt tcccaacgat caagggcagt
 3001 tacatgatcc cccatgttgt gcaaaaaagc ggtagctcc ttcggctctc cgatcgttgt
 3061 cagaagtaag ttggccgcag tgttatcact catggttatg gcagcactgc ataattctct
 3121 tactgtcatg ccatccgtaa gatgcttttc tgtgactggg gagtactcaa ccaagtcatt
 3181 ctgagaatag tgtatgcggc gaccgagttg ctcttgcccg gcgtcaatac gggataatac
 3241 cgcgccacat agcagaactt taaaagtgtc catcattgga aaacgttctt cggggcgaaa
 3301 actctcaagg atcttaccgc tgttgagatc cagttcgatg taaccactc gtgcacccaa
 3361 ctgatcttca gcacctttta ctttcaccag cgtttctggg tgagcaaaaa caggaaggca
 3421 aatgcccga aaaaagggaa taagggcgac acggaaatgt tgaatact **ca t**actcttct
 3481 ttttcaatat tattgaagca tttatcaggg ttattgtctc atgagcggat acatatttga
 3541 atgtatttag aaaaataaac aaataggggt tccgcgcaca tttccccgaa aagtgccacc
 3601 tgacgtctaa gaaaccatta ttatcatgac attaacctat aaaaataggc gtatcacgag
 3661 gccctttcgt ctgcgcggtt tcggtgatga cggtgaaaac ctctgacaca tgcagctccc
 3721 ggagacggtc acagcttgct tgtaagcggg tgccggggagc agacaagccc gtcagggcgc
 3781 gtcagcgggt gttggcgggt gtcggggctg gcttaactat gcggcatcag agcagattgt
 3841 actgagagtg caccatatgg acatattgtc gttagaacgc ggctacaatt aatacataac
 3901 cttatgtatc atacacatac gatttaggtg acactata

Beta-lactamase