

## Recombinant Gal4-VP16 activator

For *in vitro* Transcription (IVT) assays with RNAPII

**NOTE 1:** For general notes on recombinant protein expression see the *E.coli* lab protocols section.

**NOTE 2:** Gal4-VP16 is a potent transcriptional activator (see: Sadowski *et al* (1988) *Nature* **335**:563). Contains the DNA-binding domain of yeast Gal4 fused to a highly acidic portion of the herpes simplex virus protein VP16. Potent activator when added to RNAPII IVT reactions in which the template (**eg.** (pGAL4CG<sup>-</sup>; **F348 / b1457**) contains a Gal4-binding site near the promoter.

See Keogh Notebook: **pH7-143, pH7-150, pH8-9.**

1. Transform expression plasmid pGAL4-VP16 (BE367 / b2911; **Amp<sup>R</sup>**) (from Pat Nakatani) into *E.coli* strain BL21/DE3 (or BL21/DE3 pLysS). This fusion protein is HIS6 tagged internally (GAL4-his6-VP16).
2. Grow 100ml culture overnight at 37°C (LB<sup>Amp</sup>; LBA). Sub-innoculate next morning 1/100 into 2L LBA and grow to OD<sub>600</sub> ≈ 0.6 at 37°C. Chill the cultures 15 minutes on ice and induce with 1mM IPTG at RT° for two hours.

**NOTE 3: Supplement culture media with 20μM ZnSO<sub>4</sub>** (aids solubility). All steps from this point are performed at 4°C, and all buffers contain 10μM ZnSO<sub>4</sub> and protease inhibitors as standard.

3. Harvest cells by centrifugation, wash pellet with ice-cold ddH<sub>2</sub>O and transfer to a 50ml Falcon tubes in (can freeze pellets -80°C at this point). Resuspend in 30mls **buffer A** (4°C) and sonicate on ice (3 x 20s pulses 50% output, 40 duty cycle). Add Triton X-100 to 1% and mix by inversion.
4. Clarify by centrifugation (10K, 15m, 4°C) and transfer the supernatant to a 50ml falcon tube on ice. Add 2ml Ni-NTA resin (Qiagen) (pre-equilibrate in **buffer A**) and incubate at 4°C overnight with gentle rolling.
5. Collect beads by centrifugation (1.5K, 5m, 4°C) and wash in batch with 20mls **buffer A** and 20mls **buffer E / 20mM Imidazole**. Collect after each resuspension by centrifugation and discard the supernatants (keep some to test by Western analysis). Transfer the beads to a disposable column and elute with **buffer E / 200mM Imidazole**. Collect 500μl fractions.
6. Analyze fractions by Biorad Dc and SDS-PAGE (**GAL4-his6-VP16 ≈ 30 kD**). Estimate the cleanliness of the prep (how many bands visible on commassie staining SDS-PAGE) and yield per L.

**7.1 Option 1:** Pool positive fractions and dialyze against 2 x 1L *buffer E / 0mM Imidazole* (2 hrs per, 4°C). Compare concentration (by Biorad Dc) before and after dialysis (to ensure no holes in the membrane) and store in small aliquots at -80°C.

**7.2 Option 2 (clean up further by FPLC: may not be necessary):** Pool positive fractions and load onto a 1.5ml Mono-Q column (Biorad: Flow Rate (FR) 0.5ml/min). Wash with 5mls *buffer E / 0mM Imidazole*. Apply a 20ml linear gradient of 100-700mM NaCl in *buffer E / 0mM Imidazole* (Gal4-VP16 elutes at ~350mM NaCl). Pool positive fractions, estimate concentration by Biorad Dc and store in aliquots at -80°C.

**NOTE 4:** It may be necessary to clean the Mono-Q column (if it hasn't been used in a while or the **back pressure looks suspiciously high**). If so reverse the flow and successively wash (FR: 0.5ml/min) with:

- 500µl 2M NaCl
- 5ml ddH<sub>2</sub>O
- 500µl 2M NaOH
- 5ml ddH<sub>2</sub>O
- 500µl 1% TFA
- 5ml ddH<sub>2</sub>O

Then reverse the flow to normal and equilibrate in *buffer E / 0mM Imidazole*.

## **Buffers:**

### ***Buffer A -***

10mM Tris.HCl (pH 8.0)  
 10% Glycerol  
 500mM NaCl  
 0.1% Tween 20  
 10mM Imidazole  
 10mM β-ME  
 10µM ZnSO<sub>4</sub>  
 Protease inhibitors as standard (incl. 1mM PMSF)

### ***Buffer E -***

20mM HEPES.KOH (pH 7.5)  
 20% Glycerol  
 100mM NaCl  
 1mM EDTA  
 1mM DTT  
 10µM ZnSO<sub>4</sub>  
 Protease inhibitors as standard (incl. 1mM PMSF)  
 Imidazole as indicated

**For a manuscript methods section:**

Recombinant GAL4-VP16 fusion protein: E. coli strain BL21/DE3 was transformed with pGAL4-VP16 (kind gift of Pat Nakatani). A 2L culture was grown at 37°C to  $OD_{600} \approx 0.5$  and induced by the addition of 1mM IPTG. Cells were further cultured for two hours at room temperature in the presence of 20 $\mu$ M ZnSO<sub>4</sub>. All steps from this point on were performed at 4°C and all buffers contained 10 $\mu$ M ZnSO<sub>4</sub> and protease inhibitors as standard. Lysates were prepared by sonication in lysis buffer (10mM Tris.Cl pH8.0, 500mM NaCl, 10% Glycerol, 10mM  $\beta$ -ME, 0.1% Tween-20, 10 $\mu$ M ZnSO<sub>4</sub>, 10mM Imidazole, 1mM PMSF), Triton X-100 added to 0.1% and extracts clarified by centrifugation. The soluble extract was incubated at 4°C with gentle rolling overnight with 2ml Ni-NTA resin (Qiagen). The resin was applied to a column, washed with 15mls buffer E (20mM HEPES-KOH pH7.6, 100mM NaCl, 20% Glycerol, 1mM EDTA, 1mM DTT, 10 $\mu$ M ZnSO<sub>4</sub>, 20mM Imidazole, 1mM PMSF) and eluted with buffer E / 200mM Imidazole. Positive fractions were pooled, applied to a 1.5ml Mono-Q column (Biorad) and washed with 5mls buffer E (no imidazole). A 20ml linear gradient of 100-700mM NaCl in buffer E was applied where GAL4-VP16 elutes at  $\approx$  350mM NaCl. Positive fractions were pooled and stored in aliquots at -80°C.