

S.cerevisiae fractionation

NOTE 1. In addition to fractionating *S.cerevisiae*, this protocol should work for purifying chromatin / nucleosomes suitable for direct analysis or affinity purification of binding factors.

Variant of protocol used in: Keogh *et al* (2005) *Cell* **123**:593-605
Keogh *et al* (2006) *Genes Dev* **20**:660-5
Keogh *et al* (2006) *Nature* **439**:497-501

1. Grow 100mls – 1L cells in appropriate medium to $OD_{\lambda 600} < 1.0$ (exponential phase, etc). All steps below are with pre-chilled solutions on ice unless stated otherwise. Volumes given are for cells harvested from ≈ 50 mls with $OD_{600} < 1.0$. For this pellet size you can perform the sucrose gradient fractionation step in a 2ml ependorff tube (this is important – with a fixed angle microfuge, a 1.5ml ependorff is not recommended).
2. Collect cells by centrifugation. Successively wash with 10mls ddH₂O and 10mls SB in a 14ml U-bottom Falcon tube. Store pellets @ -80°C .
3. Thaw pellets on ice. Successively wash with 1.5mls each PSB and SB in a 2ml ependorff tube (flash spin to 9K in a benchtop microfuge).
4. Resuspend completely in 1ml SB. Add 125 μ l Zymolase 20T (10mg/ml in SB). Remove 10 μ l and add to 1ml 1% SDS; determine $OD_{\lambda 600}$, blanking with 10 μ l SB in 1ml 1% SDS. Incubate RT $^{\circ}$ with rotation. Determine $OD_{\lambda 600}$ every 30min, continuing incubation until it drops > 80% (should take about one hour). NOTE: Some mutants may be more zymolase-resistant than WT: if digestion is suspiciously slow at one hour add another 75 μ l Zymolase 20T and continue incubation.
5. When digestion is complete add 1ml SB, collect by centrifugation (microfuge: 2K, 5m, 4°C) and repeat SB wash (**NB. Be careful – spheroblasts are more fragile than cell-walled yeast**).
6. Resuspend spheroblasts gently with 500 μ l EBX (0.25% Triton X-100 + Phenol red). Accounting for the pellet volume, add Triton X-100 to 0.5% final to lyse outer cell membrane. Place on ice for 10min with gentle mixing every few minutes (* **aliquot / westerns; TOTAL**). Layer over 1ml NIB. Centrifuge 12K, 15m, 4°C .
7. Should now be three layers : (i) cloudy reddish (because of the phenol red) upper (cytoplasmic fraction - * **aliquot / westerns; CYTOPLASM**), (ii) possibly a white lipid (?) interphase (it may be on the side of the tube), and (iii) nuclear pellet: glassy, white, and the stuff you want.
8. Remove the supernatant and gently resuspend nuclei in 500 μ l EBX (0.25% Triton X-100). Accounting for the pellet volume, add Triton X-100 to 1% final to lyse nuclear membrane. Place on ice for 10min with gentle mixing every few minutes (* **aliquot / westerns; NUCLEI**).

9. Collect chromatin by centrifugation (15K, 10m, 4°C). Remove S/N (* **aliquot / westerns; NUCLEOPLASM**). There are numerous options at this point depending on what you want to do.

(i) **Salt stability:** Wash the chromatin from above (with EBX) and split into numerous aliquots in 1.5ml endpdroff tubes. Wash individual preps twice with EBX –100 to –500 (all mM NaCl; EBX is usually 100mM; In my hands most histones are kicked out between –250 and –500, so choose steps accordingly). Resuspend pellet in 50µl Tris pH8, and add 100µl 2x SDS-PAGE loading b°. Incubate 95°C for 5 minutes, microfuge (15K, 5m, RT°C) and collect supernatant to a new 1.5ml endpdroff tube. Run 10-15µl samples on SDS-PAGE gels to see if your modification of interest is relatively sensitive.

(ii) **Solubilize nucleosomes for immunoprecipitation:** Wash the chromatin from above (with EBX) and resuspend to 500µl with EBX. Sonicate 3 x 20s as per ChIP output (12) with one min on ice between pulses (this was determined empirically: MCK notebook pH13.152). Collect insoluble by centrifugation (15K / 10m / 4°C). PCI extract an aliquot of the supernatant and run on an EtBr gel to see how much DNA was released and average size. Also run samples on SDS-PAGE gel to see how many histones were released. Determine concentration by Biorad Dc to balance different samples and bind onto beads for further experiments.

Other useful reagents = Control Abs to confirm efficient fractionation

1. Rabbit α-Rpt1 (19S proteasome component), use 1/1000; mainly cytoplasmic
2. Rabbit α-Htz1 (Upstate 07-718), use 1/400, nuclear
3. Rabbit α-H3 (Abcam ab1791), use 1/750; Nuclear
3. Rabbit α-eIF-5a (gift from R.S. Zitomer; SUNY Albany), use 1/4000; cytoplasmic

Buffers

SB 1M Sorbitol
20mM Tris.Cl pH 7.4

PSB 20mM Tris.Cl pH 7.4
2mM EDTA
100mM NaCl
10mM β-ME

EBX 20mM Tris.Cl pH 7.4
100mM NaCl
0.25% Triton X-100
15mM β-ME

Protease inhibitors as STD

Phosphatase inhibitors if desired (1M NaF = 1000x; 0.5M Na₃VO₄ = 1000x, 1M β-Glycerophosphate = 200X)

0.005% Phenol red (from 0.5% stock (1000x) in PBS) as indicated. Note that in addition to marking the upper layer, this is a pH indicator (yellow 6.8 – 8.2 red).

NIB 20mM Tris.Cl pH 7.4
 100mM NaCl
 1.2M Sucrose
 15mM β -ME
 +/- 50mM Na-butyrate
 Protease inhibitors as STD
 Phosphatase inhibitors if desired

(If you decide to try MNase to release the nucleosomes rather than Sonication)

1x MNase b° 20mM Tris.Cl pH 8.0
 100mM NaCl
 0.1% Triton X-100
 5mM MgCl₂
 2mM CaCl
 10mM β -ME
 Protease inhibitors as STD
 Micrococcal Nuclease @ 1 – 50U as determined empirically.

2xSDS-PAGE b° 120mM Tris pH 6.8
 2% SDS
 20% Glycerol
 0.2% bromophenol blue
 200mM DTT