

S.cerevisiae transformation

NOTE 1: This is one of those protocols where the numbers are a little flexible – and it should really be idiot proof : if it didn't work you did something really stupid.

NOTE 2: Are you sure that the strain to be transformed doesn't already contain the resistance marker or auxophic marker you're about to cover ?

1. Grow 10ml culture to $OD_{600} \approx 1.0$ (assuming WT-like growth, 10AM: inoculate a loopful of a relatively fresh strain into YPD, should be right by 4PM). Collect by centrifugation: 5min, 2000rpm in a tabletop centrifuge (should get $\approx 200\mu\text{l}$ pcv) (see **NOTE 3**: this step may not be required).

NOTE 3: For everyday transformations (plasmid into strain), you can use a small loopful (20 μl pcv; packed cell volume) of *S.cerevisiae* (direct from a plate less than a month old and stored appropriately at 4°C) and 2 μl of Qiagen grade DNA (miniprep from a pRS series plasmid: 1.5mls overnight, 50 μl elution $\approx 100\text{ng/ml}$). This should give you ≥ 50 colonies. If you want targeted homologous recombination you need to step it up a bit (as in **Step #1** above).

2. Discard supernatant. Collect pellet by resuspending in 750 μl ddH₂O. Transfer to a sterile ependorff on ice. Collect by centrifugation in a table-top microfuge (a flash-spin up to 10K should be enough). Discard the supernatant and resuspend pellet in 1ml LiTE. Collect by flash-spin centrifugation. Resuspend in $< 500\mu\text{l}$ LiTE (enough for about 10 transformations).
3. Prepare ependorff tubes containing the DNA of interest (5 μl PCR product, excised fragment, whatever). Add 5 μl sheared, salmon sperm DNA (10mg/ml; this is usually stored at -20°C in 500 μl aliquots. Boil an aliquot (5m, 100°C) and store at 4°C for repeated use). Add 50 μl of the cell suspension from **Step #2** and mix well. Add 250 μl PLATE. Mix well.
4. Incubate 30°C for 15m with occasional agitation (flick the tube if you remember). Heat shock 42°C for 60-90m with occasional agitation (try and flick the tube once or twice).
5. Flash spin to collect cell pellet. Note this will take a little longer because of the PEG. Also the quality of the pellet will be different, with a significant amount along the side of the tube. Remove the supernatant and resuspend in 250 μl ddH₂O. Plate all as appropriate.

NOTE 4: If transforming a resistance marker do not plate directly. Rather resuspend pellet from **Step #5** in 500 μl YPD. Incubate (30°C, 3-4 hours) to allow gene expression and plate all as appropriate. If using KanMX (G418 selection; 200 $\mu\text{g/ml}$) it is recommended to replica the plate after 24 hours (it reduces the background significantly). With Nat (Noursethrin; 100 $\mu\text{g/ml}$) this is not necessary.

Buffers:

1xLiTE 0.1 M Lithium Acetate (1M stock)
 1x TE (10 mM Tris-HCl pH7.9, 1 mM EDTA)

PLATE 40% Polyethylene Glycol 3500 or less (50% Stock)
 0.1 M Lithium Acetate (1M Stock)
 1x TE (10 mM Tris-HCl pH7.9, 1 mM EDTA)